

SimplyP Animal pathogens DNA/RNA Extraction Kit



Product Name SimplyP Animal pathogens DNA/RNA Extraction Kit

Packing Size 24 Tests/box; 48 Tests/box; 96 Tests/box

Usage

Extract viral DNA from tissue homogenate supernatant and serum, plasma, whole blood and ascites samples.

Principle

The liquid sample is processed by the Lysis Buffer to break the virus and denature most of the protein to release the viral nucleic acid. At the same time, special polymer materials are used to selectively adsorb the viral nucleic acid, and then the viral nucleic acid can be obtained through washing and elution operations.

Kit Components

Cat#	BSC70T1	BSC70S1	BSC70M1	Components
Kit Size	24 Tests	48 Tests	96 Tests	
Lysis Buffer	7.2 mL	14.4 mL	28.8 mL	Salt and Tris Buffer
Wash Buffer	※3 mL	※6 mL	※12 mL	Low-salt solution
Elution Buffer	10 mL	10 mL	10 mL	DNase-free H ₂ O
Spin Columns	24	48	96	Plastic parts and nucleic acid adsorption film
Collection Tube	24	48	96	Plastic parts
Handbook	1	1	1	/

Notes : Please add 12mL, 24mL, 48mL of absolute ethanol to the ※ 3mL, ※ 6mL, ※ 12mL Wash Buffer before use.

Reagents to be supplied by user

Absolute ethanol (AR)

Storage and transportation

For research use only



SimplyP Animal pathogens DNA/RNA Extraction Kit



- 1. The kit can be transported at room temperature.
- 2. The kit should be stored at room temperature
- 3. The kit has demonstrated stability of 12 months when stored at room temperature.

Equipment to be supplied by user

- 1. Microcentrifuge capable of 14,000rpm;
- 2. Metal bath or Water bath;
- 3. Vortex mixer

Sample Requirements

If the volume of the liquid sample is less than 100μL, you can add an appropriate volume of PBS buffer or saline to make the total volume reach 100μL.

Protocol

- 1. Sample Preparation

Sample	Sample Preparation
Whole Blood	Take 300μL of whole blood directly, or centrifuge at 3000g for 10 minutes, then take 100μL of plasma for sample extraction.
Plasma, Serum	Take 100μL of plasma or serum sample directly for sample extraction.
Saliva	Centrifuge at 12,000g for 2 minutes, and take 100μL supernatant for sample extraction.
Feed	Add an appropriate amount of saline or PBS to the feed, grind it thoroughly, and take 100 μL of supernatant for sample extraction after centrifugation.
Disease Tissue	Take an appropriate amount of tissue sample, mix it with saline or PBS in a ratio of 1:10, grind into a homogenate, centrifuge at 12,000g for 2 minutes, and take 100μL of supernatant for sample extraction.
Swab	Add 300-400μL of normal saline or PBS to the swab collection tube, vortex the tube, and then take 100μL of supernatant for sample extraction.

For research use only



2. Reagent Preparation

Store the kit at room temperature. If the room temperature is too low, it may cause salt out in the Lysis Buffer. If this happens, place the Lysis Buffer in a 56°C water bath for 10 minutes to ensure that the salt crystal in the solution is fully dissolved.

3. Operation of Sample Extraction

- a) Take 100µL of the processed sample (please fill up to 100µL with PBS buffer or saline for samples less than 100µL), add 300µL of Lysis Buffer, shake vigorously for 30 seconds, and incubate at 80°C for 2 minutes.

Note: For samples that are difficult to lyse, the warm bath time can be extended to 5 minutes.

- b) Transfer the above mixture to the purification column, centrifuge at 12,000 rpm for 1 minute, and discard the collection tube.
- c) Transfer the purification column to a new collection tube, add 400µL of Wash Buffer to the purification column, centrifuge at 12,000 rpm for 2 minutes, and discard the collection tube.
- d) Transfer the purification column to a new 1.5mL centrifuge tube, add 50µL of Elution Buffer to the purification column, and incubate at room temperature for 1 minute.
- e) Centrifuge at 12,000 rpm for 1 minute, and discard the purification column. The liquid in the 1.5 mL centrifuge tube contains DNA.

Note: If the extracted DNA is not used immediately, please store it at -20°C or -80°C for later use.

For research use only

SimplyP Animal pathogens DNA/RNA Extraction Kit



Explanation of test results

This kit is suitable for plasma, serum, ascites and other liquid samples; tissue samples require high-sensitivity PCR detection reagents.

Limitations of the test method

Sample size: extract serum and plasma no more than 100μL;

Sensitivity: a high-sensitivity PCR detection reagent is required.

Performance Indicators

The extracted nucleic acid is tested by the high sensitivity HBV DNA detection reagents that the sensitivity of the kit reaches 100 IU/mL. The linear range reaches 100 IU/mL ~ 10⁷ IU/mL. The result is repeatedly tested and confirmed by national standard quality-controlled product

Notes

1. Before use, please check whether there are crystals in the Lysis Buffer. If the solution is crystallized because of the low temperature, please place the Lysis Buffer in a 56℃ warm bath. During the incubation, you can keep shaking to dissolve the crystals before use.
2. Add absolute ethanol to the Wash Buffer according to the label.

Basic Information

Hangzhou Bioer Technology Co.,Ltd

Address: No.1192 Bin'An Rd, Binjiang District, Hangzhou, Zhejiang Province, China

Tel: 0571-87774567 Fax: 0571-87774553

Web: www.bioer.com.cn

Zip Code: 310053

Aftersales Service Provider: Hangzhou Bioer Technology Co.,Ltd

For research use only