



N-BIOTEK
Leading Biotechnology

Adipose derived Stem Cell Isolation Kit



STEMPIA

FDA

CE


**ISO
9001**

**ISO
14001**


What is SVF?

Stromal Vascular Fraction (SVF) is a component of the lipoaspirate obtained from liposuction of excess adipose tissue. Lipoaspirate, the waste product of liposuction (cosmetic surgery), contains a large population of stem cells called **adipose derived stem cells (ASCs)**, which share a number of similarities with bone marrow stem cells, including the capacity for multilineage differentiation.

Autologous SVFs containing ASCs are currently being used in clinical settings for various orthopedic applications for human patients without any serious side effects.




SVF Treatment Process



STEMPIA *Easy and Effective SVF Isolation system!*

STEMPIA is the trade mark of NBIOTEK's SVF isolation kit. This kit consist of two specially designed tubes for isolation and concentration of SVF derived from adipose tissue. STEMPIA offer easier and shorter process for SVF isolation using less consumables and tools as compared with typical manual SVF isolation process using pipets and conical tubes. Save both cost and time for SVF isolation with STEMPIA.

SVF Isolation Procedure with STEMPIA



Total Process time with STEMPIA™ is less than 1hour

Required Consumable

Tools, Consumables for STEMPIA



VS


Tools, Consumables for MANUAL



Components

1st tube for Isolation


2nd tube for washing




Specification

	1 st Tube	2 nd Tube
Dimension	140*38mm	
Capacity	Max.80mL (Max Adipose tissue 40mL)	Max. 50mL
Tube Material	Tritan, Polypropylene, Medical Silicon *Tritan : environmental material where endocrine-disrupting bisphenol-A(BPA) not found	
Packing	Gamma Sterilized and sealed packing	
Certificate	FDA(Class I), CE(Class I) , ISO9001, ISO14001	

Performance Data



► Avg. 94.1% Viability from 6 Test Groups



► Avg. SVF Volume from 31 test groups
1.33 x 10⁶ cells in 1ml Adipose Tissue

“High Viability & Yield rate despite simpler & shorter process”

Equipment for SVF process



LED Cell Activator



Centrifuge

Multi Work Station
<Centrifuge + Shaker in Clean Bench>



N-BIOTEK
Leading Biotechnology

Website : www.n-biotek.com
E-mail : export@n-biotek.com Phone : +82 32 321 2100
402-803 Techno-Park, Pyeongcheon-ro 655,
Wonmi-Gu Bucheon-Si, GyeongGi-Do, KOREA