



Bio Basic Inc.

A world leader in serving science

CERTIFICATE OF ANALYSIS

Product	Water, sterile
Product Code	WW1002
Grade	Ultra Pure
CAS#	N/A
Mol. Formula	H ₂ O
Mol. Weight	18.02
Lot No	

Test Items	Specifications	Actual Results
Appearance	Clear Solution	
Resistivity (25 °C)	10 -18.2 MΩ·cm	
pH (at 25°C)	Check & Report	
RNase, Dnase and Protease	Not Detected	

Storage: 18-25°C.

OF 21 Rev 1

V-2018.1

20 Konrad Cres. Markham Ontario L3R 8T4 Canada
Tel: (905) 474 4493, (800) 313 7224 Fax: (905) 474 5794
Email: order@biobasic.com Web: www.biobasic.com