

10xTBE Buffer
A00269



Product Description: 10XTBE, 1000 mL

Total Product Size: 1000ml

Individual Container Size: 1000ml

Number of Containers: 1

Shipping Conditions: RT

UNSPSC Code: 41105323UNSPSC

Category:
Electrophoresis Buffers

Appearance: Clear solution

Conductivity (20%, Water) @25oC: Report

pH (20%, Water) @25oC: 8.2-8.4

Note: A single-strength (1X) solution contains 0.089M Trisbase, 0.089M borate and 0.002M EDTA.

Storage: RT. Protect from moisture.