

**Product name:**

PBS (Phosphate buffered saline),

**Catalog Number:** PD0435

**Application:**

Phosphate Buffered Saline tablets have been specifically developed for use in immunoassay procedures.

**Reconstitution:**

One tablet dissolved in 100 ml of deionized water yields 1X PBS buffer. The final 1X solution contains: 2.7 mM potassium chloride, 137 mM sodium chloride and 10 mM Phosphate Buffer, pH: 7.3-7.5 at 25 °C (1 Tablet in 100ml water).

**Storage Temperature:** Room temperature.