

TBE Buffer (Tris-Borate-EDTA)
10x Premixed Powder

A0024



Product Description:

TBE buffer (Tris-Borate-EDTA), Premix powder:

Total Product Size: A0014: 1PK (1L)
A0024: 1PK (4L)

Number of Containers: 1

Shipping Conditions: RT

UNSPSC Code: 41105323UNSPSC

Category: Electrophoresis Buffers

Storage: Store at room temperature. Protect from moisture.

Appearance: White powder

pH(1X Concentration): 8.1-8.5

Solubility (1X Concentration): PASS

Note: One pack makes 4liters of 10X TBE buffer or 40 liter 1X TBE buffer.