

Link LABEL™

BlueTooth Enabled Labeler



LABeler™ Lab Printer

The new LinkLabel™ Lab Printer brings a new level of freedom to laboratory labeling. MTC's free, downloadable LabelPro™ App allows users to endlessly customize and specialize their labels. Use the app to easily change text, font, positioning, size, and label length. The LinkLabel also has the ability to print batches, serial numbers, and a variety of barcodes.

With 12 available label cassettes, the customization doesn't end on your cell phone. Choose from multiple colors of tape, including transparent. Three label heights are available for different sized vessels. 6mm tapes are great for tissue culture plates, 12mm tapes are ideal from microtubes, and the extra tall 24mm tapes are ideal for large tubes and flasks. Special "ultra-low" label tapes are also available for storage in -80°C environments or for liquid nitrogen applications.



Label cassettes

- › Compatible with BlueTooth devices
- › Free app allows for limitless customization
- › Hundreds of preloaded scientific symbols
- › Labels in white, clear, colored and Ultra Low Temperature
- › Prints barcodes, QR cods, serial numbers and batches
- › Compact – 5.25 x 2.75 x 6.75in. (LxWxH)

+
INFO

Solicítala a
dir.comercial@gentechbio.com

MTC Bio

www.gentechbio.com · contacto@gentechbio.com

[f](#) gentechbio [t](#) @gentechbio



GENTECH
Genetics & Technology