

# TUBOS DE MICROCENTRÍFUGA CON TAPÓN DE ROSCA SURESEAL™

MTC Bio

Los microtubos estériles con tapón de rosca SureSeal™ se encuentran fabricados de polipropileno ultra claro, estos microtubos con tapón de rosca incluyen tapas ensambladas con juntas tóricas de silicona que aseguran un sello a prueba de fugas de -90 ° C a + 121 ° C. Todos los microtubos con tapón de rosca SureSeal™ son estériles, esterilizables en autoclave, libres de ARNasa, DNasa y pirógenos. Los microtubos y las tapas que los acompañan son de polipropileno y se pueden centrifugar hasta 22.000 xg.

- Polipropileno ultra claro
- Esterilizado con rayos gamma
- A prueba de fugas
- Estéril, RNasa, DNasa y de pirógenos

Cat NO.	Descripción del producto	Cantidad
C3150-S	ScrewSeal™ 1.5mL (12.8 x 45mm), sterile, w/ O-ring, conical bottom, non-graduated, caps assembled, 100 per bag, 10 bags per case	1000
C3150-SL	ScrewSeal™ 1.5mL (13 x 47mm), w/ loop caps, sterile, w/ O-ring, conical bottom, non-graduated, caps assembled, 100 per bag, 10 bags per case	1000
C2230	ScrewSeal™ 2.0mL (12.8 x 46mm), w/ screw-cap, w/ O-ring, molded graduations, self-standing, caps assembled, 100 per bag, 10 bags per case	1000
C2231-W	ScrewSeal™ 2.0mL (12.8 x 46mm), w/ screw-cap, w/ O-ring, molded graduations, self-standing, caps assembled, individually wrapped	500
C2230-AM	ScrewSeal™ 2.0mL (12.8 x 47.4mm), w/ amber screw-cap, w/ O-ring, molded graduations, self-standing, amber caps assembled, 100 per bag, 10 bags per case	1000
C3220-SG	ClearSeal™ 2.0mL (13.1 x 46.8mm), sterile, w/ O-ring & attached caps, printed graduations, self-standing, 50 per sterile resealable bag, 20 bags per case	1000
C3220-SGL	ClearSeal™ 2.0mL (13 x 49.2mm), sterile, w/ O-ring & attached loop-caps, printed graduations, self-standing, 50 per sterile resealable bag, 20 bags per case	1000