

# TUBOS DE MICROCENTRÍFUGA ESTÉRILES SURESEAL S

MitcBio

Los tubos SureSeal S™ se fabrican en condiciones de sala limpia, utilizando solo el copolímero BioClear™, un polipropileno avanzado que maximiza la claridad, la fuerza centrífuga (probado hasta 28.000 xg) y el rango de temperatura utilizable (-90 ° C a + 121 ° C). La tapa de tubo de gran tamaño (con superficie de marcado mate), combinada con el exclusivo cierre de anillo y ranura, mejora la seguridad del sello hermético a los líquidos, al tiempo que minimiza las fuerzas de bloqueo necesarias para una apertura / cierre más ergonómico de las tapas de tubo.

- Todos los tubos SureSeal S están libres de ARNasa, DNasa y pirógenos.
- Estéril por radiación de haz de electrones
- Bolsa Self-Standing – con fondo reforzado

Cat NO.	Descripción del producto	Cantidad
C2000	Microtube w/ cap, 1.5ml, clear, sterile, w/ self-standing bag & Stop-Pops™	500
C2000-AST	Microtube w/ cap, 1.5ml, assorted colors (B, R, G, O, P, Y), sterile w/ self-standing bag & Stop-Pops™	500
C2000-50	Microtube w/ cap, 1.5ml, clear, sterile, packed in bags of 50 tubes, 10 bags per unit	500
C2000-50AST	Microtube w/ cap, 1.5ml, assorted colors (B, R, G, P, Y), sterile, packed in bags of 50 tubes, 10 bags per unit (2 bags of each color)	500
C2000-B	Microtube w/ cap, 1.5ml, blue, sterile, w/ self-standing bag & Stop-Pops™	500
C2000-R	Microtube w/ cap, 1.5ml, red, sterile, w/ self-standing bag & Stop-Pops™	500
C2000-G	Microtube w/ cap, 1.5ml, green, sterile, w/ self-standing & Stop-Pops™	500
C2000-O	Microtube w/ cap, 1.5ml, orange, sterile, w/ self-standing bag & Stop-Pops™	500
C2000-P	Microtube w/ cap, 1.5ml, purple, sterile, w/ self-standing bag & Stop-Pops™	500
C2000-Y	Microtube w/ cap, 1.5ml, yellow, sterile, w/ self-standing bag & Stop-Pops™	500
C2000-OB	Microtube w/ cap, 1.5ml, opaque black, sterile, w/ self-standing bag & Stop-Pops™	500
C2002	Microtube w/ cap, 2.0ml, clear, sterile, w/ self-standing bag & Stop-Pops™	500
C2007	Microtube w/ cap, clear, sterile, 0.5ml, w/ self-standing bag	500
C2007-AST	Microtube w/ cap, 0.5ml, assorted (B, R, G, O, P, Y), sterile, w/ self-standing bag	500
C2086	Stop-Pop™ Locking Clips, w/ breakaway lifting tabs for 1.5mL tubes	100
C4010	PolyPestle™ double-sided pestle, for use with cell strainers or microcentrifuge tubes, sterile, individually wrapped	50