

Cetyltrimethylammonium Bromide CTAB



CERTIFICATE OF ANALYSIS

Product	CETYLTRIMETHYLAMMONIUM BROMIDE (CTAB)
Grade	High Purity Grade
Product Code	CB0108
Formula	C ₁₉ H ₄₂ NBr
MW	364.09
CAS#	57-09-0
LOT#	

Test Items	Specifications	Results
Appearance	White to almost white powder	
Color (Photoelectric colorimetry)	≤100 Hazen	
Moisture (K.F.)	≤0.8%	
Assay (Titration)	≥99.0%	

Storage: Store at room temperature. Protect from moisture