

TBE Buffer Tris Borate EDTA



CERTIFICATE OF ANALYSIS

| | |
|--------------|---|
| Product | TBE BUFFER(TRIS-BORATE-EDTA)10X PREMIXED POWDER |
| Grade | BIOTECH GRADE |
| Product Code | A0024 |
| Description | A single strength (1X) solution contains 0.089M Tris base, 0.089M Borate & 0.002M EDTA. |
| Lot No | |

| Test Items | Specifications | Results |
|-------------------------------|----------------|---------|
| Appearance | White powder | |
| pH(1X Concentration) | 8.1-8.5 | |
| Solubility (1X Concentration) | PASS | |

Note: One pack makes 4liters of 10X TBE buffer or 40 liter 1X TBE buffer.

Storage: Store at room temperature. Protect from moisture.

QF 21 Rev 2015