

PCR Strip Tubes CP0110



GSBIO®

无锡国盛生物工程有限公司

FEATURES

- High clarity thin film for PCR / qPCR
- Good sealing performance
- Easy to peel off, low residue
- Free from DNase, RNase and non pyrogenic

Specification

Cat NO.	Color	Tube+Cap Heigt(mm)	Product Description	PCS/ PACK	PACK/ CASE	PCS/ CASE
CP100	Clear	20.7+0.8	0.2ml Single PCR tubes	1000	10	10000
CP0100	Clear	15.7+0.8	0.1ml 8-Strip Tubes& Optical Flat Caps	125	10	1250
CP0101	White	15.7+0.8	0.1ml 8-Strip Tubes& Optical Flat Caps	125	10	1250
CP0110	Clear	20.7+0.8	0.2ml 8-Strip Tubes& Optical Flat Caps	125	10	1250
CP0111	White	20.7+0.8	0.2ml 8-Strip Tubes& Optical Flat Caps	125	10	1250
CP4000	Clear	15.7+0.8	0.1ml 12-Strip Tubes& Optical Flat Caps	125	10	1250
CP4001	White	15.7+0.8	0.1ml 12-Strip Tubes& Optical Flat Caps	125	10	1250
CP4010	Clear	20.7+0.8	0.2ml 12-Strip Tubes& Optical Flat Caps	125	10	1250
CP4011	White	20.7+0.8	0.2ml 12-Strip Tubes& Optical Flat Caps	125	10	1250