

Pipette Controller

Pipette NEO



- ⊠ User Friendly pipette controller, 167g only
- ⊠ Long-life Lithium Ion battery is installed. About 1700 times dispensing are possible by full charge (4 hours).
- ⊠ Dial-typed adjustable speed controller for your choice.
- ⊠ UV resistant for clean bench work.
- ⊠ Self- standing type without holder.
- ⊠ Fully autoclavable spare parts. (Nosecone, Sylicon Holder, and Filter)
- ⊠ Complies with CE and RoHS.



Specifications

Cat. No	Description
00-PMNEO	Pipette Mate NEO

Batteries	Lithium Ion(3.6V, 600mAh)
Charger	100/240V 50/60Hz,0.2A
Filter	20mmφ 0.2μm (Hydrophobic)
Pipette Types	plastic or glass from 1 to 100mL *Plastics and glass pipettes not sold by Nichiryo

Cat. No	Description
00-PM-N003	Battery cover
00-PM-N005	Head fitting
00-PM-N013	Tube holder
00-PM-N100	Filter, 0.2μm
00-PM-N024	Travel charger 5V US
00-PM-N030	Travel charger base unit
00-PM-N025	Travel charger UK adaptor*
00-PM-N026	Travel charger AU adaptor*
00-PM-N027	Travel charger EU

* UK and AU adaptors need a travel charger base unit.